
BD SmartSite™ Bag

REF 200100-04 REF 200250-04 REF 200500-04

- en** SmartSite™ Bag
- fr** Poche SmartSite™
- es** Bolsa SmartSite™
- it** Sacca SmartSite™
- de** SmartSite™-Beutel
- sv** SmartSite™ påse
- nl** SmartSite™-zak
- bg** Плик SmartSite™
- da** SmartSite™-pose
- fi** SmartSite™-pussi
- pt** Saco SmartSite™
- cs** Vak SmartSite™
- el** Ασκόζ SmartSite™
- et** Kott SmartSite™
- hu** SmartSite™ tasak
- hr** Vrećica SmartSite™
- lv** SmartSite™ maisis
- pl** Worek SmartSite™
- ro** Pungă SmartSite™
- ga** Mála SmartSite™
- sk** Vrecko SmartSite™
- sl** Vrečka SmartSite™
- it** SmartSite™ maišelis
- mt** Borża SmartSite™

Gilero, LLC
4319 S Alston Ave, Suite 100
Durham, NC 27713 USA

MedEnvoy Global BV
Prinses Margrietplantsoen 33
Suite 123
2595 AM The Hague
The Netherlands

Distributed by
BD Switzerland Sàrl, Route de Crassier 17,
Business Park Terre-Bonne, Batiment A4, 1262 Eysins Switzerland

BD, the BD Logo, and SmartSite™ are trademarks of Becton, Dickenson and Company or its affiliates. © 2019 BD. All rights reserved.

INDICATIONS FOR USE
The BD SmartSite™ Bag is an empty container used for administration of intravenous solutions to the patient using an intravascular administration set. Medication transfer in and out of the container is done using aseptic technique.

- PRECAUTIONS**
- DO NOT use if the outer package is damaged.
 - DO NOT use a needle to inject solution into the SmartSite™ Bag. Using a needle will damage the SmartSite™ Needle-Free Valve, causing it to leak.
 - DO NOT fill bag above nominal volume.
 - Protect SmartSite™ Bag from freezing.
 - DO NOT leave slip luer connectors unattended.
 - Use solution within 24 hours of preparation or according to storage requirements in drug package insert or facility protocol, whichever is shorter.
 - The SmartSite™ Bag is only intended to be used by trained medical professionals.
 - For customers in the European Union (EU), notify BD and your regional Competent Authority in the event of a serious incident.
 - Do not reuse. Product reuse can be unsafe because of the risk of: cross-infection, endotoxin reaction, patient injury due to material failure and/or chemical reaction due to cleaning or sterilization agents

FILLING DIRECTIONS – USE ASEPTIC TECHNIQUE

- Remove SmartSite™ Bag from overwrap packaging by tearing along perforation line.
- Prior to every access, swab the top of the SmartSite™ Bag Needle-Free Valve with 70% isopropyl alcohol (2-5 seconds) and allow to dry (approximately 30 seconds). Dry time depends on temperature, humidity, and ventilation of the area.
- Attach an ISO 80369-7 male luer lock syringe to the SmartSite™ Bag Needle-Free Valve and inject solution. Repeat, as needed, to add additional fluid to the SmartSite™ Bag. Dispose of syringe per facility protocol.
- Mix SmartSite™ Bag contents per facility protocols.

Note: The contents of the SmartSite™ Bag may be administered through either the SmartSite™ Needle-Free Valve or the spike port. Proceed to step 5a if administering through the spike port, or step 5b if administering through the SmartSite™ Needle-Free Valve.



SPIKE PORT DELIVERY DIRECTIONS – USE ASEPTIC TECHNIQUE

- Twist the spike port cap to remove, and dispose per facility protocol.
- Insert IV set spike into exposed pierceable septum until a secure connection is achieved.
- Administer contents of SmartSite™ Bag per facility protocol.
- Once contents of the SmartSite™ Bag have been administered, dispose of SmartSite™ Bag per facility protocol.

SmartSite™ DELIVERY DIRECTIONS – USE ASEPTIC TECHNIQUE

- Prior to every access, swab the top of the SmartSite™ Bag Needle-Free Valve with 70% isopropyl alcohol (2-5 seconds) and allow to dry (approximately 30 seconds). Dry time depends on temperature, humidity, and ventilation of the area.
- Attach an ISO 80369-7 male luer lock connector to the SmartSite™ Bag Needle-Free Valve.
- Administer contents of SmartSite™ Bag per facility protocol.
- Once contents of the SmartSite™ Bag have been administered, dispose of SmartSite™ Bag per facility protocol.

PVC, DEHP or Natural Rubber Latex are not part of the material formulation.

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